

S/S BOX 35 X 18 X 8 CM

Réf. AI35158



S/S BOX 35 X 18 X 8 CM



Réf. AI35158

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019659
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	18,0
Length (cm) :	35,0	Length packaging (cm) :	35,0
Height (cm) :	8,0	Height packaging (cm) :	8,0
Weight (kg) :	1,250	Weight packaging (kg) :	1,250
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:49 - This article is available on our website : www.gsh-med.fr with the following reference AI35158

All characteristics are given for information only and may be changed