

S/S BOX PERFORATED 40 X 20 X 9 CM PERFORATED TOP / BOTTOM

Réf. AI40209P



S/S BOX PERFORATED 40 X 20 X 9 CM PERFORATED TOP / BOTTOM

Réf. AI40209P

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040067889
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	21,5
Length (cm) :		Length packaging (cm) :	43,0
Height (cm) :		Height packaging (cm) :	11,0
Weight (kg) :	1,400	Weight packaging (kg) :	1,510
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:49 - This article is available on our website : www.gsh-med.fr with the following reference AI40209P