


THERMOSCAN PRO 6000 TYMPANIC THERMOMETER (SMALL CRADLE)

Réf. DWPRO6K01

WelchAllyn®



CE
0297



THERMOSCAN PRO 6000 TYMPANIC THERMOMETER (SMALL CRADLE)

WelchAlllyn®

Réf. DWPRO6K01

Infrared, tympanic thermometer with disposable probe covers

- Measurement time: 2–3 seconds
- ExacTemp™ technology helps support reliability and accuracy by detecting the stability of probe placement during measurement
- Displayed temperature range: 20.0–42.2°C (68–108°F)
- Operating temperature: 10.0–40.0°C (50–104°F)
- Accuracy: ±0.2°C (±0.4°F) for the range 35.0–42.0°C (95–107.6°F)
- Display type: Backlight liquid crystal display, four digits plus special icons
- Automatic power off: Approximately 10 seconds after last measurement
- Cradle options: Small (holds 20 covers) or large (holds 40 covers)
- Other features: Timer, memory recall and C/F conversion buttons
- Power supply: 2 × 1.5 V type AA (Alkaline) batteries or Rechargeable Battery Pack
- Dimensions: 152 mm x 44 mm x 33 mm (6" x 1,7" x 1,3")

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0297	Barcode	00732094309003
Guarantee (month):	36	Customs Code	902519
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :	3,3	Width packaging (cm):	16,0
Length (cm) :	4,4	Length packaging (cm) :	25,0
Height (cm) :	15,2	Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,500
SPECIFIC STORAGE CONDITION			
Temperature min	-25	Temperature max	55
Humidity min	15	Humidity max	95



30/03/2022 23:40 - This article is available on our website : www.gsh-med.fr with the following reference DWPRO6K01

All characteristics are given for information only and may be changed

2/2