

PERSONAL WEIGHING SCALE KERN MPD_M - 250 KGS - GRAD 200 G - CLASS III

Réf. EKMPD_M01_CIII

KERN



CE
0297



PERSONAL WEIGHING SCALE KERN MPD_M - 250 KGS - GRAD 200 G - CLASS III



Réf. EKMPD_M01_CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297	Barcode	3661040187402
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	49,0	Width packaging (cm):	
Length (cm) :	36,5	Length packaging (cm) :	
Height (cm) :	12,0	Height packaging (cm) :	
Weight (kg) :	10,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 00:05 - This article is available on our website : www.gsh-med.fr with the following reference

All characteristics are given for information only and may be changed

PERSONAL WEIGHING SCALE KERN MPD_M - 250 KGS - GRAD 200 G - CLASS III



Réf. EKMPD_M01_CIII

COMPOSITION

Code	Descriptive	Quantity
EKCERTI06	EU-Verification (CE) electr.balances 50 kg - 350 kg	1
EKMPD_M01	Personal floor scale MPD	1