

## COLUMN TO ELEVATE THE DISPLAY DEVICE(950 MM)FOR KERN BOBP,MPT,MWS,UOB

***KERN***

Réf. EKMWS02



## COLUMN TO ELEVATE THE DISPLAY DEVICE(950 MM)FOR KERN BOBP,MPT,MWS,UOB



Réf. EKMWS02

FEATURES			
Medical device	No		
Marquage		Barcode	4045761096714
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :	9,5	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMWS02