

WEIGHING PLATFORM MWS 300-400 KG 100 G - CLASS III MEDICAL

Réf. EKMWS03-CIII

KERN



CE
0297



WEIGHING PLATFORM MWS 300-400 KG 100 G - CLASS III MEDICAL



Réf. EKMWS03-CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297	Barcode	3661040063300
Guarantee (month):	36	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	106,0	Width packaging (cm):	
Length (cm) :	125,5	Length packaging (cm) :	
Height (cm) :	6,9	Height packaging (cm) :	
Weight (kg) :	28,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:06 - This article is available on our website : www.gsh-med.fr with the following reference EKMWS03-CIII

WEIGHING PLATFORM MWS 300-400 KG 100 G - CLASS III MEDICAL



Réf. EKMWS03-CIII

COMPOSITION

Code	Descriptive	Quantity
EKMWS03	Wheelchair platform scale MWS	1
EKCERTI12	EU-Verification (CE) electr.balances >350 kg-1500 kg	1
EKTRANSPORT02	Freight Cost Big Products	1