

REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 5 MM), 15CM

Réf. IA08415



REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 5 MM), 15CM



Réf. IA08415

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040006437
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,1	Width packaging (cm):	
Length (cm) :	16,6	Length packaging (cm) :	
Height (cm) :	0,4	Height packaging (cm) :	
Weight (kg) :	0,026	Weight packaging (kg) :	0,026
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/11/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference IA08415

All characteristics are given for information only and may be changed