

REDON INTRODUCING NEEDLE CURVED (DIAM. 3 MM), 15CM

Réf. IA08509



REDON INTRODUCING NEEDLE CURVED (DIAM. 3 MM), 15CM



Réf. IA08509

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040006451
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,5	Width packaging (cm):	
Length (cm) :	15,5	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,018	Weight packaging (kg) :	0,018
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/05/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference IA08509

All characteristics are given for information only and may be changed