

GUYON CATHETERS WITHOUT GUIDE 15 PCS (N°6 TO N°20)

Réf. IB02121



GUYON CATHETERS WITHOUT GUIDE 15 PCS (N°6 TO N°20)



Réf. IB02121

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040039275
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	
Length (cm) :	26,5	Length packaging (cm) :	
Height (cm) :	1,5	Height packaging (cm) :	
Weight (kg) :	0,010	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:46 - This article is available on our website : www.gsh-med.fr with the following reference IB02121

All characteristics are given for information only and may be changed