

OPERATING SCISSORS B/B STRAIGHT 14 CM, FOR LEFT-HANDED

Réf. IC11414G



OPERATING SCISSORS B/B STRAIGHT 14 CM, FOR LEFT-HANDED



Réf. IC11414G

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe Ir | |
| Marquage | CE | Barcode | 3661040083308 |
| Guarantee (month): | 12 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 5,4 | Width packaging (cm): | |
| Length (cm) : | 14,2 | Length packaging (cm) : | |
| Height (cm) : | 0,5 | Height packaging (cm) : | |
| Weight (kg) : | 0,050 | Weight packaging (kg) : | 0,050 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 15:59 - This article is available on our website : www.gsh-med.fr with the following reference IC11414G

All characteristics are given for information only and may be changed