

## B/B SCISSORS 14 CM STRAIGHT COLLIN LOCK

Réf. IC11614



# B/B SCISSORS 14 CM STRAIGHT COLLIN LOCK



Réf. IC11614

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040009889
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,4	Width packaging (cm):	
Length (cm) :	14,2	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,051	Weight packaging (kg) :	0,051
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:00 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IC11614

All characteristics are given for information only and may be changed