

## JESCO UNIVERSAL SCISSORS BLACK 17 CM

Réf. IC14400



# JESCO UNIVERSAL SCISSORS BLACK 17 CM



**Réf. IC14400**

Jesco (or Gesco) scissors.  
 Put them in every first aid kit!  
 Widely used in surgical procedures and their presence is necessary in first aid kits at school, on work sites, in company pharmacies or in retirement homes.  
 Stainless steel cutting blades.  
 Universal scissors allowing to cut easily.

Available in different sizes : yellow, black, green, blue, red  
 Universal size : 17 cm at the shortest, 19 cm at the longest  
 Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.  
 The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.  
 Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.  
 100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040010052
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,8	Width packaging (cm):	20,0
Length (cm) :	17,0	Length packaging (cm) :	13,0
Height (cm) :	6,2	Height packaging (cm) :	
Weight (kg) :	0,063	Weight packaging (kg) :	0,080
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



20/06/2022 21:05 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IC14400

All characteristics are given for information only and may be changed