

JARS POLI 1L WITH LID FOR ASPIRET, ASKIR 20,30

Réf. BCASPI78

CA-MI



CE
0124



JARS POLI 1L WITH LID FOR ASPIRET, ASKIR 20,30

CA-MI

Réf. BCASPI78

FEATURES

Medical device No

Marquage CE0124

Barcode 3661040151083

Guarantee (month): 24

Customs Code 90330090

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) : 12,0

Width packaging (cm):

Length (cm) : 21,0

Length packaging (cm) :

Height (cm) : 12,0

Height packaging (cm) :

Weight (kg) : 0,200

Weight packaging (kg) : 0,200

SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 0

Humidity max 85



30/03/2022 22:54 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI78

2/2