

DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800

WelchAllyn®

Réf. DWPRO6K02



CE

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
I



DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800



Réf. **DWPRO6K02**

This covers probe helps protect your ThermoScan Pro 6000 ear thermometer. Single use for better efficiency and impeccable hygiene, this protective cap is perfectly suited to the Braun thermometer. Its shape is studied to suit children and babies. Colour: transparent

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	00732094321043
Guarantee (month):	0	Customs Code	9025900090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	30,0
Height (cm) :		Height packaging (cm) :	8,0
Weight (kg) :		Weight packaging (kg) :	0,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:40 - This article is available on our website : www.gsh-med.fr with the following reference DWPRO6K02

All characteristics are given for information only and may be changed