

INTERFACE CABLE FOR PC FOR MODELS KERN ABT, CGB

KERN

Réf. EKMP510



INTERFACE CABLE FOR PC FOR MODELS KERN ABT, CGB



Réf. **EKMPS10**

FEATURES			
Medical device	No		
Marquage		Barcode	4045761014824
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	2,7
Length (cm) :		Length packaging (cm) :	6,0
Height (cm) :		Height packaging (cm) :	3,5
Weight (kg) :		Weight packaging (kg) :	80,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:06 - This article is available on our website : www.gsh-med.fr with the following reference EKMPS10

All characteristics are given for information only and may be changed