

VACUUM GAUGE FOR ASKIR

Réf. BCASPI150



VACUUM GAUGE FOR ASKIR



Réf. BCASPI150

FEATURES			
Medical device	No		
Marquage		Barcode	0000000000949
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,5
Length (cm) :		Length packaging (cm) :	5,5
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,042
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI150