

POWER CABLE FOR NEW ASKIR 30

Réf. BCASPI170



POWER CABLE FOR NEW ASKIR 30



Réf. BCASPI170

FEATURES			
Medical device	No		
Marquage	CE0124	Barcode	3661040031552
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,158
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	5	Humidity max	85



31/03/2022 00:30 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI170