

**CONNICAL CONNECTOR Ø 8-9-10  
MM-ASPIRET,ASKIR12-20-30-230/12BVR,36BR,  
E**



Réf. BCASPI168



## CONNICAL CONNECTOR Ø 8-9-10 MM-ASPIRET,ASKIR12-20-30-230/12BVR,36BR, E

CA-MI

Réf. BCASPI168

### FEATURES

Medical device No

Marquage

Barcode

3661040160832

Guarantee (month): 0

Customs Code

Sterile

No

### INFORMATIONS

Main unit

Packaged unit

Width (cm) : 7,5

Width packaging (cm):

Length (cm) : 9,0

Length packaging (cm) :

Height (cm) : 1,5

Height packaging (cm) :

Weight (kg) : 0,100

Weight packaging (kg) : 0,010

### SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 0

Humidity max 85



15/11/2022 12:44 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI168

2/2