

CARDIOMATE PC ECG SPENGLER



Réf. 260039



CARDIOMATE PC ECG SPENGLER



Réf. 260039

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446042059
Guarantee (month):	24	Customs Code	9018110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,6	Width packaging (cm):	
Length (cm) :	6,7	Length packaging (cm) :	
Height (cm) :	2,4	Height packaging (cm) :	
Weight (kg) :	0,130	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/05/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference 260039

All characteristics are given for information only and may be changed