

## CARDIOMATE PC ECG BLUETOOTH SPENGLER

Réf. 260040



# CARDIOMATE PC ECG BLUETOOTH SPENGLER



Réf. **260040**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446042066
Guarantee (month):	24	Customs Code	9018110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,7	Width packaging (cm):	
Length (cm) :	12,6	Length packaging (cm) :	
Height (cm) :	2,4	Height packaging (cm) :	
Weight (kg) :	0,130	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/05/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 260040

All characteristics are given for information only and may be changed