

OXYFROG - DIGITAL FINGERTIP PULSE OXIMETER SPENGLER

Réf. 221056



OXYFROG - DIGITAL FINGERTIP PULSE OXIMETER SPENGLER



Réf. 221056

Ideal for measuring the pulse of children

High measurement accuracy:

- High sensitivity sensor
- Suitable for children from 10 kg to 40 kg

Easy to use and convenient:

- Compact & lightweight
- Automatic shut-off after 5 seconds
- Easy to use: 1 button

Additional information:

- Comes with a wrist strap and 2 AAA LR03 batteries

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0123	Barcode	06945040105442
Guarantee (month):	24	Customs Code	90181910
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,7	Width packaging (cm):	7,3
Length (cm) :	5,9	Length packaging (cm) :	11,7
Height (cm) :	3,5	Height packaging (cm) :	4,9
Weight (kg) :	0,135	Weight packaging (kg) :	0,135
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/05/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference 221056

All characteristics are given for information only and may be changed