

SPENGLER TB100 GLUCOMETER COMPLETE KIT

Réf. 745000



SPENGLER TB100 GLUCOMETER COMPLETE KIT



Réf. 745000

FEATURES			
Medical device	Yes	DMDIV (diagnostic In Vitro)	
Marquage	CE0537	Barcode	3700446041649
Guarantee (month):	24	Customs Code	9027808000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,8	Width packaging (cm):	12,0
Length (cm) :	9,2	Length packaging (cm) :	16,5
Height (cm) :	1,9	Height packaging (cm) :	5,0
Weight (kg) :	0,600	Weight packaging (kg) :	0,220
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/05/2022 12:35 - This article is available on our website : www.gsh-med.fr with the following reference 745000

All characteristics are given for information only and may be changed