

Réf. 745010



Réf. 745010

FEATURES			
Medical device	Yes	DMDIV (diagnostic In Vitro)	
Marquage	CE0537	Barcode	3700446041656
Guarantee (month):	24	Customs Code	90183190
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	11,5
Length (cm) :		Length packaging (cm) :	16,4
Height (cm) :		Height packaging (cm) :	5,2
Weight (kg) :		Weight packaging (kg) :	0,220
SPECIFIC STORAGE CONDITION			
Temperature min	4	Temperature max	30
Humidity min	0	Humidity max	0



31/03/2022 03:28 - This article is available on our website : www.gsh-med.fr with the following reference 745010

All characteristics are given for information only and may be changed