

SOFTWARE PC-ECG FOR ECG CARDIOMATE (WITH PHYSICAL KEY)



Réf. 260038

SOFTWARE PC-ECG FOR ECG CARDIOMATE (WITH PHYSICAL KEY)



Réf. **260038**

| FEATURES | | | |
|----------------------------|-----------|-------------------------|-------|
| Medical device | No | | |
| Marquage | Barcode | | |
| Guarantee (month): | 0 | Customs Code | |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 03:27 - This article is available on our website : www.gsh-med.fr with the following reference 260038

All characteristics are given for information only and may be changed