

LIAN METAL / GREY COTTON CUFF / ADULT (M)

Réf. 521318



CE
0459



LIAN METAL / GREY COTTON CUFF / ADULT (M)



Réf. 521318

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446008840
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,750
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



06/09/2022 12:40 - This article is available on our website : www.gsh-med.fr with the following reference 521318

All characteristics are given for information only and may be changed