

## 5X PROTECTIVE WORKING COVER MPE-A01 FOR MPE AND MTA

**KERN**

Réf. EKMPE04



## 5X PROTECTIVE WORKING COVER MPE-A01 FOR MPE AND MTA



Réf. **EKMPE04**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040173559
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	5,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMPE04

All characteristics are given for information only and may be changed