

## KERN MTA WEIGH PLATFORM, RESTRAINT BAR-MEDICAL CLASS III

Réf. EKMTA01\_CIII

**KERN**



# KERN MTA WEIGH PLATFORM, RESTRAINT BAR-MEDICAL CLASS III



Réf. EKMTA01\_CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297		Barcode 3661040173580
Guarantee (month):	36		Customs Code
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	68,0	Width packaging (cm):	
Length (cm) :	78,0	Length packaging (cm) :	
Height (cm) :	110,0	Height packaging (cm) :	
Weight (kg) :	40,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:32 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMTA01\_CIII

## KERN MTA WEIGH PLATFORM, RESTRAINT BAR-MEDICAL CLASS III



Réf. EKMTA01\_CIII

### COMPOSITION

Code	Descriptive	Quantity
EKMTA01	Handrail scale MTA	1
EKCERTI12	EU-Verification (CE) electr.balances >350 kg-1500 kg	1
EKTRANSPORT02	Freight Cost Big Products	1