

FOLDING SIDEBAR AND DISPLAY HOLDER

KERN

Réf. EKMWA04



FOLDING SIDEBAR AND DISPLAY HOLDER



Réf. **EKMWA04**

FEATURES			
Medical device	No		
Marquage	Barcode	3661040173672	
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	40,0	Width packaging (cm):	
Length (cm) :	56,0	Length packaging (cm) :	
Height (cm) :	101,3	Height packaging (cm) :	
Weight (kg) :	4,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:32 - This article is available on our website : www.gsh-med.fr with the following reference EKMWA04

All characteristics are given for information only and may be changed