

XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III

Réf. EKMWS07_CIII

KERN



CE
0297



XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III



Réf. EKMWS07_CIII

FEATURES			
Medical device	Yes	Classe IIa	Accuracy measurement
Marquage	CE0297		Classe III
Guarantee (month):	36	Barcode	3661040173740
Sterile	No	Customs Code	
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	86,0	Width packaging (cm):	
Length (cm) :	150,0	Length packaging (cm) :	
Height (cm) :	6,8	Height packaging (cm) :	
Weight (kg) :	40,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:32 - This article is available on our website : www.gsh-med.fr with the following reference EKMWS07_CIII

All characteristics are given for information only and may be changed

XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III



Réf. EKMWS07_CIII

COMPOSITION

Code	Descriptive	Quantity
EKMWS07	MWS-L	1
EKTRANSPORT02	Freight Cost Big Products	1
EKCERTI06	EU-Verification (CE) electr.balances 50 kg - 350 kg	1