

## MATRIX NEEDLE PRINTER PRINTOUT SUITABLE FOR A LONGTERM STORAGE

**KERN**

Réf. EKACCE06



## MATRIX NEEDLE PRINTER PRINTOUT SUITABLE FOR A LONGTERM STORAGE



Réf. **EKACCE06**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040067063
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:01 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKACCE06