

CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM

Réf. DUCRYO03



CE
0344



CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM



Réf. **DUCRYO03**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0344	Barcode	8717484005835
Guarantee (month):	0	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	7,5
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	14,5
Weight (kg) :		Weight packaging (kg) :	0,300
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	25
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 00:00 - This article is available on our website : www.gsh-med.fr with the following reference DUCRYO03