

PAPER ROLLS FOR PRINTER 911-013 (10 PIECES) WIDTH 57 MM

KERN

Réf. EKACCE09



PAPER ROLLS FOR PRINTER 911-013 (10 PIECES) WIDTH 57 MM



Réf. **EKACCE09**

FEATURES			
Medical device	No		
Marquage	Barcode	3661040077901	
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	1,000	Weight packaging (kg) :	1,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:01 - This article is available on our website : www.gsh-med.fr with the following reference EKACCE09

All characteristics are given for information only and may be changed