

POWER SUPPLY (CH,EURO,UK) 12V, 500MA, INPUT 100V-240V (BEFORE 2017)

KERN

Réf. EKACCE11



POWER SUPPLY (CH,EURO,UK) 12V, 500MA, INPUT 100V-240V (BEFORE 2017)



Réf. **EKACCE11**

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040126807
Guarantee (month):	36	Customs Code	85044090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,400	Weight packaging (kg) :	0,400
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:01 - This article is available on our website : www.gsh-med.fr with the following reference EKACCE11

All characteristics are given for information only and may be changed