

STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Réf. AD25126



STERILIZING STAINLESS STEEL BOX 25X12X6 CM



Réf. AD25126

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040013770
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,0	Width packaging (cm):	14,8
Length (cm) :	25,0	Length packaging (cm) :	27,0
Height (cm) :	6,0	Height packaging (cm) :	8,0
Weight (kg) :	0,726	Weight packaging (kg) :	0,744
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:42 - This article is available on our website : www.gsh-med.fr with the following reference AD25126

All characteristics are given for information only and may be changed