

BAR UMBILICAL CLAMP 8 CM

Réf. IC30108



BAR UMBILICAL CLAMP 8 CM



Réf. IC30108

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040010922
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,5	Width packaging (cm):	
Length (cm) :	8,0	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,016	Weight packaging (kg) :	0,016
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:07 - This article is available on our website : www.gsh-med.fr with the following reference IC30108

All characteristics are given for information only and may be changed