

STEVENS SCISSORS STRAIGHT 18 CM

Réf. IC03418



STEVENS SCISSORS STRAIGHT 18 CM



Réf. IC03418

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040083070
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,2	Width packaging (cm):	5,2
Length (cm) :	20,3	Length packaging (cm) :	20,3
Height (cm) :	0,5	Height packaging (cm) :	0,5
Weight (kg) :	0,059	Weight packaging (kg) :	0,059
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:52 - This article is available on our website : www.gsh-med.fr with the following reference IC03418

All characteristics are given for information only and may be changed