

STITCH SCISSORS STRAIGHT 16 CM

Réf. IC03616



STITCH SCISSORS STRAIGHT 16 CM



Réf. IC03616

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040083124
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,2	Width packaging (cm):	
Length (cm) :	14,9	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,040	Weight packaging (kg) :	0,042
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:52 - This article is available on our website : www.gsh-med.fr with the following reference IC03616

All characteristics are given for information only and may be changed