


THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

Réf. DWPRO6K08

WelchAllyn®



CE
0297



THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

Réf. **DWPRO6K08**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0297	Barcode	3661041009826
Guarantee (month):	36	Customs Code	902519
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,3	Width packaging (cm):	16,0
Length (cm) :	4,4	Length packaging (cm) :	25,0
Height (cm) :	15,2	Height packaging (cm) :	15,0
Weight (kg) :	0,440	Weight packaging (kg) :	0,460
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 00:10 - This article is available on our website : www.gsh-med.fr with the following reference DWPRO6K08

THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

Réf. DWPRO6K08

COMPOSITION

Code	Descriptive	Quantity
DWPRO6K07	Charging Base Station for ThermoScan Pro 6000 Tympanic Thermometer	1
DWPRO6K01	ThermoScan Pro 6000 tympanic thermometer (small cradle)	1