


PERSONAL SCALES WITH COLUMN FOR CONVENIENT, ERGONOMICALLY-OPTIMISED (MAX. 250KG)

Réf. EKMPK01

KERN



CE
0297



PERSONAL SCALES WITH COLUMN FOR CONVENIENT, ERGONOMICALLY-OPTIMISED (MAX. 250KG)



Réf. EKMPK01

Ideal measuring instrument for elderly care, paramedics, fitness clubs or in the private sector and can handle weighing overweight persons weighing up to 250 kg

- The display device is easy to operate through the column which is firmly fixed to the scale and the weight can be read easily, even under poor lighting conditions
- Sturdy glass weighing plate for secure, comfortable standing position ; Easy and hygienic cleaning
- Large high-contrast display that is easy to read (height of digits 28 mm)
- Particularly flat design facilitates access
- Easy to operate: simply tap on it – stand on it – read the weight
- Integrated AUTO-OFF function to preserve the batteries
- Total dimensions (including column) L x W x H 350 x 420 x 955 mm
- Batteries included, 4 x 1.5V AA, service life up to 20 h
- Net weight approx. 3.8 kg

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | No | | |
| Marquage | CE0297 | Barcode | 3661040195414 |
| Guarantee (month): | 24 | Customs Code | |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 42,0 | Width packaging (cm): | 14,0 |
| Length (cm) : | 35,0 | Length packaging (cm) : | 102,0 |
| Height (cm) : | 95,5 | Height packaging (cm) : | 45,0 |
| Weight (kg) : | | Weight packaging (kg) : | 6,350 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 00:53 - This article is available on our website : www.gsh-med.fr with the following reference EKMPK01

All characteristics are given for information only and may be changed