


CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH

Réf. 570150



CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH



Réf. 570150

Suitable for all head sizes thanks to a flexible metal headband.
Indirect lighting by reflecting mirror.
Suited to lighting required in daily practice.
Built-in battery compartment, two CR 123 batteries
Supplied with carrying bag

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446011390
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,647
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference 570150

All characteristics are given for information only and may be changed