


DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER



Réf. L20747



DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER



Réf. **L20747**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446031268
Guarantee (month):	6	Customs Code	90330000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,6
Length (cm) :		Length packaging (cm) :	11,7
Height (cm) :		Height packaging (cm) :	11,4
Weight (kg) :		Weight packaging (kg) :	0,140
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:24 - This article is available on our website : www.gsh-med.fr with the following reference L20747

All characteristics are given for information only and may be changed