

TENSONIC WRIST - WRIST BLOOD PRESSURE MONITOR



Réf. 527522



CE
0459



TENSONIC WRIST - WRIST BLOOD PRESSURE MONITOR



Réf. **527522**

For regular self-monitoring. Wrist circumference 13,5-21 cm (5,3-8,5 inch). Back-lit display. 2 x 60 memories with date and time. 3 average available : all measurements of a user / last seven daystest results for morning (5:00 am ~ 8:59 am) / last seven days test results for afternoon (18:00 pm ~ 19:59 pm). Supplied with 2 LR03 AAA batteries and a storage case.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446022617
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	110,0
Length (cm) :		Length packaging (cm) :	165,0
Height (cm) :		Height packaging (cm) :	65,0
Weight (kg) :		Weight packaging (kg) :	0,288
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min	15	Humidity max	90



08/12/2022 12:37 - This article is available on our website : www.gsh-med.fr with the following reference 527522

All characteristics are given for information only and may be changed