

AUTOTENSIO BPM WIRST CUFF, BLUE SPENGLER

Réf. SPG340-



AUTOTENSIO BPM WIRST CUFF, BLUE SPENGLER



Réf. SPG340-

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446022518
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	10,0
Length (cm) :		Length packaging (cm) :	9,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	0,270
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min	15	Humidity max	90



08/12/2022 12:39 - This article is available on our website : www.gsh-med.fr with the following reference SPG340-