

## AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS

Réf. SPG340P-



# AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS



Réf. **SPG340P-**

Display of 6 AUTOTENSIO® wrist model. For occasional self-monitoring. Very easy-to-use. - Quick and quiet inflation method- 120 memories with date and time- LCD screen- Fully automatic- Purple cuff- Wrist circumference 13,5 - 21 cm (5,3-8,5 inch). Supplied with 2 LR03 AAA batteries and a storage case.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446023027
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,700
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



08/07/2022 21:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SPG340P-

All characteristics are given for information only and may be changed

## AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS



**Réf. SPG340P-**

### COMPOSITION

Code	Descriptive	Quantity
SPG340-	Autotensio BPM wirst cuff, Blue Spengler	6