

DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO

Réf. SP3100P



DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO



Réf. SP3100P

Display of 6 thermometers Tempo Duo

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446023133
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,243
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:32 - This article is available on our website : www.gsh-med.fr with the following reference SP3100P

All characteristics are given for information only and may be changed

DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO



Réf. SP3100P

COMPOSITION

Code	Descriptive	Quantity
SP3101/	Thermometre-Tempo Duo Blue	6