

## CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS

Réf. 260003



# CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS



Réf. 260003

ECG Cardiomate 3® and 6® are robust and reliable ECG. Featuring a 7" color screen tilt they allow viewing 12 lead simultaneously. Simple to use, recording and printing ECG are triggered with only one button. The Cardiomate 3® & 6® are Cardiomate equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer. The Cardiomate 3® & Cardiomate 6® are equipped with automatic interpretation and pacemaker detection. The internal battery life of 3 hours, their lightness and carrying handle make Cardiomate 3® & Cardiomate 6® the perfect companion in the office or traveling.

- Carrying handheld.
- Tilting screen : 7-inch color LCD display (resolution: 800x480 px)
- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM
- Autonomy :The operating life is 1080 ECG.
- Automatic interpretation and intervals calculation.
- Pacemaker detection.
- Network connection, USB & printer.
- Printing and rapid analysis.
- Acquisitions and display 12 simultaneous tracks.
- Alphanumeric keypad.

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE0123    | Barcode                 | 3700446025472 |
| Guarantee (month):         | 24        | Customs Code            | 9018110000    |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 24,4      | Width packaging (cm):   |               |
| Length (cm) :              | 31,0      | Length packaging (cm) : |               |
| Height (cm) :              | 6,5       | Height packaging (cm) : |               |
| Weight (kg) :              | 2,400     | Weight packaging (kg) : | 4,680         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -20       | Temperature max         | 50            |
| Humidity min               | 25        | Humidity max            | 85            |



08/12/2022 12:38 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 260003

All characteristics are given for information only and may be changed