

PURPLE SPENGLER (ML) CUFF (22-42CM) FOR AUTOTENSIO



Réf. 527343



PURPLE SPENGLER (ML) CUFF (22-42CM) FOR AUTOTENSIO



Réf. **527343**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446007522
Guarantee (month):	6	Customs Code	90262040
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	16,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,124
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference 527343