

SPO2 EXTENSION CABLE FOR MASTER PALM 2 & 3

Réf. 222237



CE
0123



SPO2 EXTENSION CABLE FOR MASTER PALM 2 & 3



Réf. 222237

SpO2 extension cable for MASTER PALM 2 & 3 oximeters with a length of 2 meters

FEATURES					
Medical device	Yes	Classe IIb	Accuracy measurement	Classe II	
Marquage	CE0123		Barcode	3700446025373	
Guarantee (month):	0		Customs Code	90330000	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		0,071	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 01:40 - This article is available on our website : www.gsh-med.fr with the following reference 222237

All characteristics are given for information only and may be changed