

SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2



Réf. 222255



SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2



Réf. **222255**

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | No | | |
| Marquage | | Barcode | 3700446030254 |
| Guarantee (month): | 12 | Customs Code | 90181990 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,084 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | 1 | Temperature max | 50 |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 01:41 - This article is available on our website : www.gsh-med.fr with the following reference 222255

All characteristics are given for information only and may be changed