


## SPO2 SENSOR ADULT FINGER CLIP FOR MASTER PALM 2

Réf. 222250



# SPO2 SENSOR ADULT FINGER CLIP FOR MASTER PALM 2



Réf. **222250**

Spo2 adult rigid sensor for MASTER PALM 2 oximeter

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0123	Barcode	3700446030209
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,103
SPECIFIC STORAGE CONDITION			
Temperature min	1	Temperature max	50
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:32 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 222250

All characteristics are given for information only and may be changed